

SW PULSATE™

Pulsation Therapy – Various modes, ease of use.



P R E C A U T I O N

Preserving

the safety

& dignity

P R E V E N T I O N

of patients

& caregivers

alike

P R O T E C T I O N

COMFORT

Automatically adjusts the mattress for changes in patient position, as well as the capability to change firmness of mattress with the touch of a button. Adjustable modes and lock-out features.

CARE

Pressure redistribution, low air loss therapy with static and pulsating modes; passive massage action aids in the increase of capillary blood flow.

CONVENIENCE

Accommodates weight capacities up to 600 lbs. Sizes to fit your needs: 35", 39" widths and step deck system options. Available in mattress or overlay, as well as pediatric sizes.



SIZEWISE

800.814.9389 | sizewise.net

SW PULSATE™

CONTROL UNIT

Liters of Air Per Minute1740
 Alternating Cycle Time N/A
 Static Low Air Loss..... Yes
 Pulsating Low Air Loss.....3, 5, or 10 minutes
 Continuous Lateral Rotation..... N/A
 Alarms..... Power Failure & Low Pressure
 Warranty 1 year
 Pressure Sensor Technology..... Yes
 Patented CPR Release..... Yes
 Dimensions.....18" W x 7 1/2" D x 11 1/2" H
 Unit Weight.....16 lbs.

MATTRESS

Weight Capacity 600 lbs.
 Mattress Size35", 39" & 48"
 Percussion Modules Add-on Feature
 Warranty 1 year
 Top Cover Material..... 200 Denier Polyurethane Coated
 Height Inflated.....7"
 Mattress Weight.....37 lbs.
 Top Cover Warranty.....90 days

ELECTRICAL

Rated Voltage 110 volts
 Rated Frequency..... 50 Hz
 Rated Input Power..... 120 volts
 Fuse Rating..... 125 volts 5 amps
 Power Failure Alarm Yes

SW PULSATE™ OVERLAY


CONTROL UNIT

Alternating Cycle Time N/A
 Static Low Air Loss..... Yes
 Pulsating Low Air Loss.....3, 5, or 10 minutes
 Continuous Lateral Rotation..... N/A
 Alarms..... Power Failure & Low Pressure
 Warranty 1 year
 Pressure Sensor Technology..... Yes
 Patented CPR Release..... Yes

MATTRESS

Weight Capacity 600 lbs.
 Mattress SizeStandard 35"
 Warranty 1 year
 Top Cover Warranty.....90 days

SW PULSATE™ CONTROL PANEL



MADE IN THE USA



800.814.9389 | sizewise.net