

# POWERDRIVE™

Freedom –  
to maneuver with ease.



## COMFORT

Adjustable height & speeds for transport.

## CARE

Cardio Chair, Trendelenburg, & Reverse Trend positions.

## CONVENIENCE

Fully electric with weight capacity up to 1,000 lbs.

P R E C A U T I O N

Preserving

the safety

& dignity

P R E V E N T I O N

of patients

& caregivers

alike

P R O T E C T I O N



SIZEWISE

800.814.9389 | [sizewise.net](http://sizewise.net)

REAL THERAPY FOR REAL PEOPLE™

## POWER DRIVE

### MECHANICAL

Battery Size.....U1 (33 amp hour gel sealed)  
Speed: Adjustable .....2.5 mph  
Range.....5-7 miles (varies with weight)  
Turning Radius .....50"  
Ground Clearance .....3"  
Max Climbing.....2"  
Tires .....foam filled  
Maximum Weight Capacity .....1000 lbs.

### ELECTRICAL

Input Voltage.....120V AC, 60Hz  
Current.....2.7 amperes  
Charger Output Voltage .....24V DC  
Classification .....Class B 1P-54

**MADE  
IN THE USA**



SIZEWISE

800.814.9389 | [sizewise.net](http://sizewise.net)