

SIZEWISE™: Pediatric Pulsation Therapy

- Dual modes of clinically effective therapy – static & pulsating

- Low friction and shear top cover

- Pulsation cycle time is adjustable from 3 to 20 minutes



REAL THERAPY for REAL PEOPLE™

- Clinically effective low air loss therapy
- Pulsating feature automatically adjusts for patient comfort and pressure redistribution
- Gentle stimulating action aids in the increase of capillary blood flow to the skin
- Top cover design reduces friction and shear
- Aids in moisture management with air flow through top cover while maintaining waterproof characteristics
- Custom designs fit any hospital or homecare crib

Clinically Effective LAL



Easy-to-use Controls



800.814.9389 | sizewise.net

SIZEWISE™: Pediatric Pulsation Therapy

PEDIATRIC PULSATE™ SPECS

CONTROL UNIT

Zoned Low Air Loss.....	Yes
Pulsating Low Air Loss.....	3, 5, or 10 minutes
Alarms.....	Power Failure & Low Pressure
Warranty.....	1 year
Pressure Sensor Technology.....	Yes
Patented CPR Release.....	Yes
Dimensions.....	18" W x 7 ½" D x 11 ½" H
Unit Weight.....	16 lbs.

MATTRESS

Maximum Weight Capacity.....	100 lbs.
Mattress Size.....	Varies
Warranty.....	2 Years Limited*
Top Cover Material.....	70 Denier Polyurethane Coated Nylon Taffeta
Height Inflated.....	4"
Mattress Weight.....	20 lbs.
Top Cover Warranty.....	90 days

ELECTRICAL

Rated Voltage.....	110 volts
Rated Frequency.....	50 Hz
Rated Input Power.....	120 volts
Fuse Rating.....	125 volts 5 amps
Power Failure Alarm.....	Yes

*Air Support limited warranty includes two (2) years on air cell welds; one (1) year on bottom cover, manifold and foam base; Ninety (90) Days – Top cover. All parts found defective within that period shall be repaired and/or replaced.



SIZEWISE

800.814.9389 | sizewise.net