

SIZEWISE™: SW Behavioral Health Bed™

Capital Purchase Only



REAL THERAPY for REAL PEOPLE™

The SW Behavioral Health Bed™ provides a user friendly design ensuring ergonomic safety for care providers in the behavioral and mental health environment. This fully functional electric medical bed can be utilized in a variety of care environments ensuring greater patient safety and comfort. The SW PowerSafe™ technology ensures electrical safety for both the care provider and patient. The SW Behavioral Health bed sets a new standard with its low deck height, the lowest bed on the market.

- Footboard nursing lock out controls accessible with universal key ensuring patient safety
- Quick disconnect hand control operates on either side of the bed for easy access and removal
- Electronic articulation for the head, knee and high/low feature
- Trendelenburg and Reverse Trendelenburg capabilities
- Contiguous molded construction of headboard, footboard and side rails
- All add-on features are bolted onto the bed frame
- Tamper-resistant encased wiring
- Ability to provide multiple therapy options on one bed frame, allowing for increased patient healing
- Integrated bed exit alarm, scale and protocol timer available
- Accommodates patients up to 850 lbs.

Universal Key Lock Out Feature



Optional Rest Secure System™



800.814.9389 | sizewise.net

SIZEWISE™: SW Behavioral Health Bed™

BED

Low Position - Deck to Floor (without Scales)	7.5"
Low Position - Deck to Floor (with Scales)	9"
High Position - Deck to Floor	30"
Sleep Deck Width	35"
Sleep Deck Length	82"
Maximum Overall Width w/ Rails	38"
Maximum Overall Length w/o Bumper	89"
Patient Weight Capacity	850 lbs

SIDE RAILS

Type	Fold Down
Included with Purchase	Yes
Bolted On	Yes

HEAD/FOOT BOARDS

Included in Purchase	Yes
Type	High Impact Plastic
Bolted On	Yes

BED ANGLES

Head Section	0° - 60°
Knee Section	0° - 50°
Max Trend	0° - 17°
Max Reverse Trend	0° - 17°


Note: Angles have 2° +/- variance

CASTER & BRAKING SYSTEM

Caster Diameter	4"
Dual Brake/Steer System	4"

ELECTRICAL SPECS

Emergency Battery back-up	Included
Safety Listing	IEC 60601
Foot End Bed Controls	Included
Power Requirements	120v, 60 cycle



Universal Key Lock Out Feature



800.814.9389 | sizewise.net

MADE IN THE USA