

# SIZEWISE™: Shuttle B Series™ Advanced



## REAL THERAPY for REAL PEOPLE™

The **Shuttle B Series™ Advanced** is the ultimate transitional item designed for use in conjunction with rehabilitation, examination, and transportation. Whether in stretcher or chair position, this fully electric product accommodates weight capacities up to 1,000 lbs. The Shuttle B Series Advanced transforms with the touch of a button from stretcher to chair position and is perfect for use with a lateral transfer device to assist with ambulating a patient. The seat assist function is complete with ambulation poles for the safety and confidence of the patient and caregiver. The seat has a 3 layer foam insert for further comfort and weight distribution and leg supports elevate independently with retractable foot rests. The Shuttle B Series Advanced also is available with a Wound Care Clinic and power drive option.

Elevating Leg Rests



Power Drive Option



Seat Assist



Tilt & Recline



Chair Up/Down



SIZEWISE

800.814.9389 | [sizewise.net](http://sizewise.net)

# SIZEWISE™: Shuttle B Series™ Advanced

## SHUTTLE

Inner Seat Width .....	31"
Seat Depth .....	19"
Seat Height from Floor .....	30" High, 21.5" Low
Stretcher Height from Floor .....	35.5" High, 27.5" Low
Overall Width .....	37.5"
Overall Length .....	69" to 75"
Degree of Tilt .....	15 degrees
Wheels .....	Independently Swiveling Casters
Product Weight .....	350 lbs
Weight Capacity .....	1,000 lbs
Electrical Rating .....	120v; 60hz; 6amps
ETL Approved .....	Yes
Accessories .....	6 Position Electronic Hand Control, Multiple Location IV Pole Holders, Storage Basket, Oxygen Tank Holder

## PATENT

U.S. Patent No. 8,028,359

## Power Drive Option: Powered by Safety™



- Safely navigate smooth floors, thick carpet and uneven terrain in chair or stretcher position.
- Propel patients (weighing up to 1,000 lbs.) up ramps, hills and door jambs with an incline of up to 8 degrees.
- Variable speeds with Auto-Off mode after 30 mins of inactivity to save on battery life.
- Center wheel design allows for easy maneuvering.

## 3 ways to MOBILIZE

### 1. Lateral Transfer



#### Featured Lateral Transfer Products:

Shuttle B Series™ Advanced  
Bari Rehab Platform™ & SW Air Transfer™

### 2. Total Lift



#### Featured Total Lift Products:

Shuttle B Series™ Advanced  
Bari Lift & Transfer™ Overhead

### 3. Ambulate



#### Featured Ambulation Products:

Shuttle B Series™ Advanced  
Bari Walker™



800.814.9389 | [sizewise.net](http://sizewise.net)