

# SIZEWISE™: SW Bedside Recliner™



## REAL THERAPY for REAL PEOPLE™

The **SW Bedside Recliner™** is designed for mobility with the bariatric patient in mind. A closed-arm style with black urethane armrests provides comfort and sophistication. Features include two rigid locking casters in front and two swivel locking casters in the rear that provide flexibility in movement while maintaining safety. The SW Bedside Recliner comes with a variety of features to support both the caregiver and patient in order to provide safety and convenience. Features like the relaxing footrest and multiple reclining positions put comfort in the hands of the patient.

- Base and handle have catalyzed finishes for added protections against alcohol and other cleaning solvents
- Strong, steel-to-steel construction provides extra long-life performance
- Standard with break away footrest & ergonomic push bar
- High density foam with moisture barrier feature on seat & back
- Quick release feature allows for removable back

### Specifications:

Overall Width	40"	Overall Depth	38.5"
Height	47"	Full Layout	73"
Seat Width	29"	Seat Height	21"
Seat Depth	21"	Product Weight	152 lbs.
Weight Capacity	750 lbs.		

MADE IN THE USA

### Extra Long-Life Construction



### Ergonomic Push Bar



SIZEWISE

800.814.9389 | [sizewise.net](http://sizewise.net)