

# SIZewise™: Evolution™



## REAL THERAPY for REAL PEOPLE™

**Introducing Evolution™** – the next generation in versatile, patient care beds. With an 850 pound capacity and class leading 9" low position, the Evolution™ is a safe option for 'at risk fall' patients. The deck height of 30.75" allows for an adjustable ergonomic working environment for care providers of all sizes and the bed can transport a patient in any high or low position. When combined with our therapeutic support surfaces, the Evolution™ can provide alternation therapy, pulsation therapy or rotational therapy with percussion and vibration to accommodate a wide array of patient needs. The Rest Secure System™ is an integrated scale, bed exit alarm and protocol timer that interfaces with the nurse call system. Caregiver lockout controls, nurse call interface, bed exit alarm, cardio chair positioning, 20 degree Trendelenberg and Reverse Trendelenberg make the Evolution™ the ideal bed for any care environment.

Cardio Chair



Trapeze Option



Rest Secure System™



Class Leading 9" Low



SIZewise

800.814.9389 | [sizewise.net](http://sizewise.net)

# SIZEWISE™: Evolution™

## available THERAPIES

### BED

|   |         |
|---|---------|
| Low Position - Deck to Floor .....      | 9"      |
| High Position - Deck to Floor .....     | 30.75"  |
| Sleep deck width .....                  | 35"     |
| Sleep deck length .....                 | 82"     |
| Maximum overall width w/ rails .....    | 38"     |
| Maximum overall length w/o bumper ..... | 89"     |
| Patient weight capacity .....           | 850 lbs |

### SIDE RAILS

|                              |                                |
|------------------------------|--------------------------------|
| Type .....                   | Fold Down                      |
| Included with purchase ..... | One Left, One Right, Head Rail |

### HEAD/FOOT BOARDS

|                            |                     |
|----------------------------|---------------------|
| Included in Purchase ..... | Yes                 |
| Type .....                 | High Impact Plastic |

### BED ANGLES

|                         |          |
|-------------------------|----------|
| Head Section .....      | 0° - 60° |
| Knee Section .....      | 0° - 50° |
| Max Trend .....         | 0° - 17° |
| Max Reverse Trend ..... | 0° - 17° |

**Note:** Angles have 2° +/- variance

### CASTER & BRAKING SYSTEM

|                               |    |
|-------------------------------|----|
| Caster Diameter .....         | 4" |
| Dual Brake/Steer System ..... | 4" |

### SCALE SYSTEM

|                             |     |
|-----------------------------|-----|
| Included in Purchase .....  | Yes |
| Weigh in any position ..... | Yes |

### ELECTRICAL SPECS

|                                 |                |
|---------------------------------|----------------|
| Emergency Battery back-up ..... | Included       |
| Safety Listing .....            | EN60601        |
| Foot End Bed Controls .....     | Included       |
| Power Requirements .....        | 120v, 60 cycle |

### MATTRESS OPTIONS

|  |          |
|--|----------|
| Low Air Loss Surface .....               | Optional |
| Alternation Therapy .....                | Optional |
| Full Body Lateral Rotation Therapy ..... | Optional |
| Pulsation Therapy .....                  | Optional |
| Non-Powered Surface .....                | Optional |
| Tri-Laminate Surface .....               | Optional |

### OPTIONS


|               |     |
|---------------|-----|
| Trapeze ..... | Yes |
|---------------|-----|

### PATENTS


U.S. Patent No. 7,013,510  
U.S. Patent No. 7,334,277




Control Panel with Independent Lockout Features




### Alternation Therapy – SW Alternate™



### Rotation Therapy – SW Rotate™



### Pulsation Therapy – SW Pulsate™



800.814.9389 | [sizewise.net](http://sizewise.net)