

# BARI-CHAIR™

**Reliable –**  
Provides a custom fit.



## COMFORT

Water resistant, breathable material.

## CARE

Extended surface calf pad and extended capacity elevating leg rests with foot plates. Swing away armrests for ease of patient transfer.

## CONVENIENCE

Seat widths ranging from 24" to 36", 750 lb. weight capacity.

P R E C A U T I O N

Preserving

the safety

& dignity

P R E V E N T I O N

of patients

& caregivers

alike

P R O T E C T I O N



**SIZEWISE**

800.814.9389 | [sizewise.net](http://sizewise.net)

### **BARI-CHAIR**

Weight Capacity .....750 lbs  
Seat Widths..... \*20", 22", 24", 26", 28", 30" (Customizable)  
Seat Depths.....\*20"-22" (Customizable)  
Seat to Floor Height..... \*20.5"-24" (Customizable)  
Back Height from Seat .....\*16" (Customizable)  
Overall Width..... 28.5", 30.5", 32.5", 34.5", 36.5", 38.5"  
Overall Height..... \*40"  
Overall Depth..... \*32"  
Folded Width..... 13" w/o Limiter  
Front Wheels.....6" x 1" Solid Rubber Casters  
Back Wheels..... \*24" Pneumatic/Spoked Wheels  
Features..... **Removable Armrests, Extended Surface Area on Calf Rests, Extended Weight Capacity on Independently Elevating Leg Rests (Swing Away or Removable), and Seat Belt**

\*=may vary

**MADE IN THE USA**



**SIZEWISE**

800.814.9389 | [sizewise.net](http://sizewise.net)