



# BARI-REHAB PLATFORM™

Adaptable –  
a bed to meet your changing needs.



## P R E C A U T I O N

Preserving

the safety

& dignity

## P R E V E N T I O N

of patients

& caregivers

alike

## P R O T E C T I O N

### COMFORT

Sizes to fit your needs: 39" x 80", 48" x 80", 39" x 86".

### CARE

Comes standard with pressure redistribution mattress; Cardio Chair, Trendelenburg, & Reverse Trend positions.

### CONVENIENCE

Fully electric with weight capacity up to 1,000 lbs. Scale, full frame trapeze, and power drive available.



**SIZEWISE**

800.814.9389 | [sizewise.net](http://sizewise.net)