

SIZEWISE™: Shuttle A Series™ Advanced



REAL THERAPY for REAL PEOPLE™

The **Shuttle A Series™ Advanced** is the ultimate transitional item designed for use in conjunction with rehabilitation, examination, and transportation. Whether in stretcher or chair position, this fully electric product accommodates weight capacities up to 650 lbs. The Shuttle A Series Advanced transforms with the touch of a button from stretcher to chair position and is perfect for use with a lateral transfer device to assist with ambulating a patient. The seat assist function is complete with ambulation poles for the safety and confidence of the patient and caregiver. The seat has a 3 layer foam insert for further comfort and weight distribution and leg supports elevate independently with retractable foot rests. The Shuttle A Series Advanced also is available with a Wound Care Clinic option and power drive option.

Elevating Leg Rests



Power Drive Option



Seat Assist

Tilt & Recline

Chair Up/Down



SIZEWISE

800.814.9389 | sizewise.net

SIZEWISE™: Shuttle A Series™ Advanced

SHUTTLE – U.S. Patent No. 8,028,359

Inner Seat Width	25"
Seat Depth.....	19"
Seat Height from Floor	30" High, 21.5" Low
Stretcher Height from Floor	35.5" High, 27.5" Low
Overall Width	32.5"
Overall Length.....	69" to 75"
Degree of Tilt.....	15 degrees
Wheels.....	Independently Swiveling Casters
Product Weight	325 lbs
Weight Capacity	650 lbs
Input Voltage.....	120V AC; 60hz; 6amps
ETL Approved	Yes
Power Drive	Optional
Accessories.....	8 Position Electronic Hand Control, Multiple Location IV Pole Holders, Holders, Oxygen Tank Holder

Power Drive Option: *Powered by Safety™*



- Safely navigate smooth floors, thick carpet and uneven terrain in chair or stretcher position.
- Propel patients (weighing up to 650 lbs.) up ramps, hills and door jambs with an incline of up to 5 degrees.
- Variable speeds with Auto-Off mode after 30 minutes of inactivity to save on battery life.
- Center wheel design allows for easy maneuvering.

MADE IN THE USA



800.814.9389 | sizewise.net